

Pressure tank 110 l

Cat. no: ZC-0110



Destiny:

Pressure tank used for water storage demineralized water.

Building:

The tank is built of a fully corrosion resistant housing made of fiberglass. The construction of the tank guarantees no contact of stored water with external environment (light, air).

- Height: 113 cm
- Diameter: 41 cm
- Total volume: 110 l
- Capacity. total at pressure 2,8 bar: 76,2 l
- Weight: 11,5 kg
- Max. pres. full collection: 3,5 bar
- Ambient temperature: 2 - 40°C
- Empty tank pressure: 0,4 - 0,5 bar