

service of the device

R

model	sediment prefilter 5µm	module A2	module H7 TOC	module H6 TOC	Capsule MF 0,2µm	UV lamp 254 nm	UV lamp 185/ 254 nm	ultrafiltration module
R 5	+	+	+	-	+	+	-	-
R5 UV	+	+	+	-	+	-	+	-
R5 UF	+	+	+	-	+	-	+	+
R10	+	+	-	+	+	+	-	-
R10 UV	+	+	-	+	+	-	+	-
R10 UF	+	+	-	+	+	-	+	+
R20	+	+	-	+	+	+	-	-
R20 UV	+	+	-	+	+	-	+	-
R20 UF	+	+	-	+	+	-	+	+
R30	+	+	-	+	+	+	-	-
R30 UV	+	+	-	+	+	-	+	-
R30 UF	+	+	-	+	+	-	+	+
Lifetime	6 months*	6 months*	2000 dm ³	5000 dm ³	12 months*	8500 h.	8500 h.	12 months*
nr kat	EO-005-10	EO-MA-12	EJ-2000-1	EJ-5000-1	EM-SP-20	EUV-254-0	EUV-185-254-0	EU-HLP-01

* If feed water parameters are different (hard water) the shell live of the spare parts might be shorter!